

**WIHS VIROLOGIC REBOUND AND RESISTANCE STUDY  
QUESTION BY QUESTION SPECIFICATIONS  
VRS 29: BLOOD SPECIMEN COLLECTION FORM**

Guidelines for completing VRS 29, "Blood Specimen Collection Form"

- A1. Record the date the blood was drawn in MM/DD/YY format.
- A2, Record the time of blood draw. Be sure to indicate AM or PM by circling the appropriate code 1= AM or 2 = PM.
- A3. Record phlebotomist's initials.
  
- B1. –B4. Listed are the test type and tubes required for VRS visits, in the order in which they should be drawn.
  - a. Indicate whether or not a tube was drawn for each test. If 'YES', go to question b. If 'NO', specify reason in the comments section. Circle '3' if a particular specimen doesn't need to be collected at this visit.
  - b. Indicate whether or not the required volume for each tube was collected. If 'YES', go to next test type, if 'NO', specify reason in comments section.